

Relax & View® Video Relaxation System for MRI

Reduced stress levels lead to:

- a relaxed and therefore co-operative patient
- in many cases no need for sedation
- faster throughput of patients
- in children the relief of boredom, with all its associated problems

Features:

- MR Conditional up to and including 3 Tesla - no artifacts on your images
- new 24" widescreen features LED backlit technology provides a bright, clear view of the displayed image with exceptional depth of colour.
- the unit is normally mounted on the wall at the rear of the scanner which alleviates the problem associated with some projection systems, where a screen is placed between patient and MRI operative.

- as the system is permanently mounted there is no setup required ensuring efficient patient throughput. There is an optional mobile stand if the system can't be mounted to the wall.
- solution for claustrophobic patients

Various studies have confirmed that our environment can greatly contribute to our well being and that a positive state of mind can greatly enhance a patients speed of healing. The Relax & View® Video Relaxation System is a unique product and helps to relieve the claustrophobia, stress and



boredom often associated with having an MRI scan. Designed and built in the UK, this system provides a complete audio / visual experience. This helps to relieve the problems associated with the MR environment and to enhance the patients feeling of well being. This also has a beneficial effect on staff. The use of advanced fibre optic technology for the delivery of the signal to the shielded screen enclosure ensures an artifact free scan under all conditions including

3T systems. The displayed image can be inverted and or reversed, so that no matter the orientation of the patient in the scanner, or the mirror configuration, the patient will always view the displayed image in its correct orientation.

The audio side can be either connected through the scanner's audio system (if fitted) or through our PureSound Audio System.

* Proven technology - 10 years in use around the world

THE RELAX & VIEW® SYSTEM IS EASILY INSTALLED IN EXISTING MRI FACILITIES
WITHOUT THE NEED TO RAMP DOWN THE MAGNET

CSP ITEM # DESCRIPTION

518-1114	Wardray MR6300, Relax & View Video Relaxation System
Options:	
518-1113	Wardray MR6003, Relax & View Stand
518-1123	One year additional warranty for Relax & View
518-1057	Wardray MR6002, MRI Vision Correction Goggles
518-1128	Wardray MR108NRR, PURESOUND Audio Relaxation System